

Mechanical Alignment of a Customized TKA Implant System as Determined by Intraoperative Computer Navigation

Gary A. Levensgood, MD

Presented at the proceedings of the 2014 British Association for Surgery of the Knee (BASK) Annual Meeting, Abstract #169

INTRODUCTION

A neutral mechanical limb alignment has been linked to increased Total Knee Arthroplasty (TKA) survivorship^{1,2}. Patient-specific cutting blocks have demonstrated close to neutral alignment on average, but have also found a varying number of outliers ($>3^\circ$) within specific brands (Otismed, Visionaire, Signature etc.) as well as between brands³⁻⁶. However, no study has investigated patient-specific cutting guides in conjunction with customized TKA implants (ConforMIS iTotal CR). Hence, the purpose of this study was to utilize computer navigation intra-operatively to measure the mechanical alignment of a patient-specific implant and instrument TKA system.

METHODS

A consecutive series of 63 patients undergoing TKA, utilizing patient-specific instruments and implants, were prospectively measured with intraoperative computer navigation. The instruments and implants were first created utilizing a pre-operative CT scan. The mechanical alignment of all patients was measured using computer navigation. The patient-specific instruments were then utilized per the manufacturer's recommendations and bone cuts were made (Figure 1). All bone cuts were recorded utilizing the navigation system as a confirmatory measurement. The patient-specific implants were then fixated and surgery completed. Final mechanical alignment was recorded with computer navigation.



Figure 1. Computer assisted navigation was used to measure the pre-op and post-op mechanical alignment intra-operatively.

RESULTS

The patient-specific instruments and implants provided neutral mechanical alignment, zero degrees, in 84.1% of patients (53/63). The remaining 10 patients each had a post-operative alignment within $\pm 2^\circ$ of neutral, with no

outliers. The average pre-operative degree of deformity for this cohort was 5.57° v. 0.19° postoperatively ($p < 0.0001$) (Table 1). The mean correction angle for this cohort was 5.70° . Additionally, no patients had extension deficits as measured with navigation post-operatively (7.5° pre-op for 40/63 patients).

Table 1: Average pre-op and post-op mechanical axis measurements

	Pre-op	Post-op
Average Deformity	$5.57^\circ \pm 6.5^\circ$	$0.19^\circ \pm 0.47^\circ$
Deformity Range	12° varus - 15° valgus	1° varus - 2° valgus
Extension Deficit	7.5° (5.58-30)	Full extension in all patients

DISCUSSION

Customized, patient-specific instruments and implants accurately restore neutral mechanical alignment as measured by intra-operative computer navigation. The patient-specific instruments aligned all patients in this cohort to within $\pm 2^\circ$ of neutral. It is well documented that the restoration of a neutral mechanical axis within $\pm 3^\circ$ is important to achieve long-term survivorship.¹

REFERENCES

1. Ritter MA, Faris PM, Keating EM, et al. Postoperative alignment of TKA. Its affect on survival. Clin Orthop Relat Res 1994;299:153
2. Green GV, Berend KR, Berend ME, et al. The effects of varus tibial alignment on proximal tibial surface strain in total knee arthroplasty; the posteromedial hot spot. J Arthroplasty 2002; 17:1033
3. Ng VY, DeClaire JH, Berend KR, et al. Improved accuracy of alignment with patientspecific positioning guides compared with manual instrumentation in TKA. ClinOrthop Relat Res 2012;470:99.
4. Nunley RM, Ellison BS, Ruh EL, et al. Are patient-specific cutting blocks costeffective for total knee arthroplasty? Clin Orthop Relat Res 2012;470:889.
5. Nunley RM, Ellison BS, Zhu J, et al. Do patient-specific guides improve coronal alignment in total knee arthroplasty? Clin Orthop Relat Res 2012;470:895.
6. Lustig S, Scholes CJ, Oussedik SI, et al. Unsatisfactory Accuracy as Determined by Computer Navigation of VISIONAIRE Patient-Specific Instrumentation for Total Knee Arthroplasty. J Arthroplasty 2013(28), 469-473.